

 Supplies **made** simple

**Interacoustics**

# ***SUPPLIES*** **HUB**

**Audiology Clinic  
Equipment Catalogue 2023**



# Contents

4-5	Disposable Ear Tips
6-7	Probe Tips
8-9	Supplies for Insert Earphones
10	Sanibel™ Infant EarCup™
11	EARTurtle™
12-13	REM Probe Tubes for HA Fitting
14-21	Wax Removal
22-24	Impression Materials & Supplies
25	Hearing Aid Fitting
26-27	Hearing Aid Maintenance & Care
28-29	Tubing
30-31	Tools
32-35	Disposable Electrodes
36-37	Electrode Cables
38-39	Patient Preparation
40-41	Infection Control
42-45	Audiometry Components
46	Balance Supplies
47	Paper
48-49	Otosopes & Otoscopy Supplies
50	Batteries & Various Items

## Abbreviation key

ABR	Auditory Brainstem Response	BAL	Balance	IMP	Impedance
AUD	Audiometry	HIF	Hearing Instrument Fitting	OAE	Otoacoustic Emissions
LATEX	Products not made with natural rubber latex	+	Medical grade silicone rubber		

Delivery charges

Order by Midday for next day delivery *(subject to availability)*

Delivery charges £10 plus VAT for all orders up to £250.

Free of charge for orders over £250

VAT

All prices quoted exclude VAT



Disposable Sanibel™ ADI Ear Tips

For use with the following acoustic immittance instruments and otoacoustic emission systems with CLEAR probe tip:  
**Amplivox:** all OtoWave versions. **Interacoustics:** Titan, AT235, AA220, Eclipse, MT10, Sera. **Maico:** EasyTymp, MI 22, MI 23, MI 24, MI 26, MI 34, MI 44, Easyscreen.



ADI Assortment Box

Part no.	Description	Price
8511445	Various sizes for Tympanometry and OAE	£43.97
8501572	Screening ear tip kit, various sizes for newborn screening	£43.97



Flanged infant ear tips

Part no.	Description	Unit	Price
8502837	Preemie (green)	Baq/25 pcs	£13.95
8502838	Preemie (green)	Baq/100 pcs	£27.90
8012963	3-5mm (red)	Baq/25 pcs	£13.95
8012964	3-5mm (red)	Baq/100 pcs	£27.90
8012965	4-7mm (blue)	Baq/25 pcs	£13.95
8012966	4-7mm (blue)	Baq/100 pcs	£27.90
8012967	5-8mm (yellow)	Baq/25 pcs	£13.95
8012968	5-8mm (yellow)	Baq/100 pcs	£27.90



Umbrella ear tips

Part no.	Description	Unit	Price
8012985	7.5mm (green)	Baq/25 pcs	£13.95
8012986	7.5mm (green)	Baq/100 pcs	£27.90
8012987	9mm (red)	Baq/25 pcs	£13.95
8012988	9mm (red)	Baq/100 pcs	£27.90
8012989	10mm (blue)	Baq/25 pcs	£13.95
8012990	10mm (blue)	Baq/100 pcs	£27.90
8012991	11mm (yellow)	Baq/25 pcs	£13.95
8012992	11mm (yellow)	Baq/100 pcs	£27.90
8012993	13mm (green)	Baq/25 pcs	£13.95
8012994	13mm (green)	Baq/100 pcs	£27.90
8012995	15mm (red)	Baq/25 pcs	£13.95
8012996	15mm (red)	Baq/100 pcs	£27.90
8012997	19mm (blue)	Baq/25 pcs	£13.95
8012998	19mm (blue)	Baq/100 pcs	£27.90
8012999	22mm (yellow)	Baq/25 pcs	£13.95
8013000	22mm (yellow)	Baq/100 pcs	£27.90



Mushroom ear tips

Part no.	Description	Unit	Price
8013001	7mm (blue)	Baq/25 pcs	£13.95
8013002	7mm (blue)	Baq/100 pcs	£27.90
8013003	8mm (yellow)	Baq/25 pcs	£13.95
8013004	8mm (yellow)	Baq/100 pcs	£27.90
8012969	9mm (green)	Baq/25 pcs	£13.95
8012970	9mm (green)	Baq/100 pcs	£27.90
8012971	10mm (red)	Baq/25 pcs	£13.95
8012972	10mm (red)	Baq/100 pcs	£27.90
8012973	11mm (blue)	Baq/25 pcs	£13.95
8012974	11mm (blue)	Baq/100 pcs	£27.90
8012975	12mm (yellow)	Baq/25 pcs	£13.95
8012976	12mm (yellow)	Baq/100 pcs	£27.90
8012977	13mm (green)	Baq/25 pcs	£13.95
8012978	13mm (green)	Baq/100 pcs	£27.90
8012979	14mm (red)	Baq/25 pcs	£13.95
8012980	14mm (red)	Baq/100 pcs	£27.90
8012981	15mm (blue)	Baq/25 pcs	£13.95
8012982	15mm (blue)	Baq/100 pcs	£27.90
8012983	19mm (yellow)	Baq/25 pcs	£13.95
8012984	19mm (yellow)	Baq/100 pcs	£27.90



DISPOSABLE SANIBEL™ AOT & ANA Ear Tips



ANA ear tips for echo-Screen®

For use with **Natus:** echo-Screen®.  
**Madsen/Otometrics:** AccuScreen® Classic.  
**Interacoustics:** Eclipse with pink probe tip

Part no.	Description	Unit	Price
8106363	4.0mm (blue)	Baq/50 pcs	£18.60
8106362	4.5mm (red)	Baq/50 pcs	£18.60
8106364	5.0mm (green)	Baq/50 pcs	£18.60
8106366	Infant tree (yellow)	Baq/50 pcs	£18.60
8106365	Mushroom (green)	Baq/50 pcs	£18.60



Ear Tips for OAE screeners

For use with **Biologic AuDx**, **Welch Allyn** OAE screener and **Natus** echo-Screen® III.



OAE ear tips, AuDx

Part no.	Description	Unit	Price
8107885	AuDx Preemie 3mm	Baq/100 pcs	£80.24
8515505	AuDx Preemie 3mm	Baq/20 pcs	£18.18
8107889	AuDx Tree 3.5-5mm	Baq/100 pcs	£80.24
8515506	AuDx Tree 3.5-5mm	Baq/20 pcs	£18.18
8107887	AuDx Infant 6mm	Baq/100 pcs	£80.24
8515507	AuDx Infant 6mm	Baq/20 pcs	£18.18
8102681	Adult OAE foam	Baq/100 pcs	£80.24
8515436	Adult OAE foam	Baq/20 pcs	£18.18
8102682	Pediatric OAE foam	Baq/100 pcs	£80.24
8515435	Pediatric OAE foam	Baq/20 pcs	£18.18
8103123	Jumbo OAE foam	Baq/100 pcs	£87.17
8515438	Jumbo OAE foam	Baq/20 pcs	£19.53



AOT ear tips

For use with **Madsen/Otometrics** AccuScreen®, Biologic Bio-logic® AuDx®, **Pathmedical:** all Sentiero products

Part no.	Description	Unit	Price
8106367	3.7mm (yellow)	Baq/50 pcs	£18.60
8106368	4.0mm (light green)	Baq/50 pcs	£18.60
8106369	4.5mm (light pink)	Baq/50 pcs	£18.60
8106370	5.0mm (light blue)	Baq/50 pcs	£18.60
8106371	4-7mm (transparent)	Baq/50 pcs	£18.60



Age	Ear tips inserted against ear canal	Ear tips inserted against entrance of ear canal	Ear tips for Echo-Screen®	Ear tips for AccuScreen®
0.5 y/o		Green pediatric	4.0mm, ANA	3.7mm, AOT
			4.5mm, ANA	4.0mm, AOT
		3-5mm, flanged	5.0mm, ANA	4.5mm, AOT
	7.5mm umbrella	4-7mm, flanged	infant tree, ANA	5.0mm, AOT
	9mm umbrella	5-8mm, flanged	mushroom, ANA	4-7mm, AOT
10 y/o	10mm umbrella	7mm mushroom	Selecting the correct ear tip	
	11mm umbrella	8mm mushroom		
	13mm umbrella	9mm mushroom		
Adults		10mm mushroom		
	15mm umbrella	11mm mushroom		
		12mm mushroom		
	19mm umbrella	13mm mushroom		
		14mm mushroom		
		15mm mushroom		
	22mm umbrella	19mm mushroom		



Probe tip for Otowave and new MT10™

Part no.	Unit	Price
8002592	1 pc	£5.66

IMP



Probe tip for Titan and easyTymp™

Part no.	Unit	Price
8001701	1 pc	£13.95

ABR

IMP

OAE



Probe tip kit

Fits Titan vers. 2 AT235™ or newer, AA222™, Sera™, Tympstar Pro™ Novus™, easyScreen™, and touchTymp™. Includes probe tip and gasket.

Part no.	Unit	Price
8501570	Set/3 pcs	£27.90

IMP



Probe tubes for ADI™ ear tips

Fits Corti™, ERO-SCAN™ and the OtoRead 11™.

Part no.	Unit	Price
8104159	Box/100 pcs	£20.97

ABR

OAE



Probe Tube Removal Tool

Fits Corti™, ERO-SCAN™ and the OtoRead 11™.

Part no.	Unit	Price
8513131	Bag/1 pc	£4.23

ABR

OAE




Probe Tubes 100pcs with Removal Tool kit

Fits Corti™, ERO-SCAN™ and the OtoRead 11™.

Part no.	Unit	Price
8513132	Bag/100 pcs + 1 pc	£23.25

ABR

OAE



Probe tip kit for OAE Ø3.3mm ear tips

Fits all Eclipse™ – use with ADI™ ear tips. Requires no re-calibration.

Part no.	Unit	Price
8012960	Bag/8 pcs	£69.75

OAE



Probe tip adaptor kit

Converts AT235™, AT235h™, AA220™, AA222™(blue/new), MI24™, MI26™, MI34, MI44™ for use with ADI™ ear tips. Includes probe tip, washer, damper. No recalibration required.

Part no.	Unit	Price
8500430	Set/3 pcs	£18.60

IMP



Remote probe tip replacements for ERO-SCAN™

Fits ERO-Scan™.

Part no.	Unit	Price
8102340	Bag/4 pcs	£27.90

OAE



Internal probe tip replacements for ERO-SCAN™

Fits ERO-SCAN™.

Part no.	Unit	Price
8102339	Bag/4 pcs	£31.62

OAE



ER10D probe tips for OtoRead™

It is recommended to clean the probe tip after every patient and to change the probe tip after 20 tests (40 ears). No calibration necessary.

Part no.	Unit	Price
8012836	Bag/4 pcs	£23.25

OAE



Internal probe tip replacements for ERO-SCAN™ Pro

Fits ERO-Scan™ Pro.

Part no.	Unit	Price
8120914	Bag/4 pcs	£31.62

OAE

Hearing Snap Probe Accessories

Accessories specifically designed to work with **SnapPROBE®**. For use with **Interacoustics Sera, Grasonstadler Novus, Maico EasyScreen.**



SnapPROBE™ Accessories box

Part no.	Unit	Price
8522461	Assortment of accessories incl. Infant tip A & C and waxguard	£44.64



Ear tip Infant A (Conical)

Part no.	Unit	Price
8522758	Bag/100 pcs	£48.19
8522759	Bag/25 pcs	£16.57



Ear tip Infant C (Narrow)

Part no.	Unit	Price
8522760	Bag/100 pcs	£48.19
8522761	Bag/25 pcs	£16.57



Probe Waxguard

Part no.	Unit	Price
8516436	1 pc	£5.58
8517889	Set 10 pcs	£51.15





### Sanibel™ audiological foam tips

Complies to all insert phones.

Part no.	Size	Unit	Price
8527912	Pediatric/ small, pink	Baq/ 50 pcs	£20.29
8527913	Pediatric/ small, pink	Baq/ 100 pcs	£27.48
8527915	Pediatric/ small, pink	Baq/ 750 pcs	£199.96
8527918	Adult/ medium, pink	Baq/ 50 pcs	£20.29
8527920	Adult/ medium, pink	Baq/ 100 pcs	£27.48
8527921	Adult/ medium, pink	Baq/ 750 pcs	£199.96
8535534	Adult/ Jumbo, pink	Baq/ 24 pcs	£15.22

ABR




### ADI ear tip adaptor

ADI ear tips compatible  
with 3A insert phones.

Part no.	Unit	Price
8500580	Baq/6 pcs	£18.60

ABR




### Sanibel™ IP30™ tube connectors

For use with IP30™ and  
E-A-RTone® 3A.

Part no.	Unit	Price
8107391	Baq/10 pcs	£14.46

ABR



### Replacement tube nipples clear for EARTone 3A

For use with sound tubes 3A.

Part no.	Unit	Price
8500140	Baq/10 pcs	£16.32

ABR




### Sanibel™ sound tubes and Sanibel™ ADI™ ear tip adaptor kit

Fits E-A-Rtone® 3A inserts and  
IP30™ insert phones. Includes  
ADI™ ear tip adaptors, 1 pc  
blue and 1 pc red 50 ohm.

Part no.	Unit	Price
8100590	Set/2 pcs	£29.76

ABR



### Yellow Tube and probe tip adaptor kit

Part no.	Unit	Price
8107396	Baq/4 pcs	£14.80

ABR



### Fastener kit for IP30™

Includes velcro neckstrap and clip.

Part no.	Unit	Price
8507410	Set	£29.59

ABR



### IP30 and EARTone sound tubes and nipples

Fit all IP30™ and E-A-RTone®  
versions. Includes tube connectors,  
1 pc blue and 1 pc red 50 ohm.

Part no.	Unit	Price
8501072	Set/2 pcs	£22.32

ABR



### Sanibel™ ear tip adaptor cleaning tool

Small brush for cleaning the Sanibel™  
probe tip adaptors before disinfection.

Part no.	Unit	Price
8500830	Baq/10 pcs	£2.37

ABR

OAE




### ABR consumable kit

Includes electrodes, insert phone  
tips, skin scrub and gaze.

Part no.	Description	Unit	Price
8109670	Adult 3A	Kit/100 tests	£325.52
8109671	Pediatric 3B	Kit/100 tests	£325.52

ABR



We know that, first and foremost, great quality is at  
the top of your list, and you'll find it in each and every  
purchase you make from the Interacoustics Supplies Hub.

Bringing you the products you need  
is what drives us every day.

If we can ever be of assistance,  
don't hesitate to let us know.

[supplieshub@interacoustics.com](mailto:supplieshub@interacoustics.com)

**Claire Hague**  
Director of Sales & Marketing



Sanibel™ Infant EarCup™  
with snap or tab electrodes

Reduce test time! Test on both ears in one screening.  
For use with a range of different providers of newborn hearing screening equipment: **Grason Stadler:** Novus. **Interacoustics:** Titan, Sera. **Maico:** Easyscreen, MB11 Classic. **Madsen/Otometrics:** Accuscreen. **Natus:** Biologic Abear, Echoscreen.

Part no.	Description	Unit	Price
8049600	Tab	Box/20 sets	£156.42
8049700	Snap	Box/20 sets	£156.42
8508848	w/o electrodes	Box/20 sets	£138.66

ABR



Sanibel Infant EarCup™  
with Hydrogel.

New Hydrogel version to ensure a more comfortable solution for the baby while testing.  
Choose any electrode type you wish to add to meet your preference.  
Order electrodes separately

Part no.	Unit	Price
8516095	Box/20 sets	£148.81



Sanibel™ Infant EarCup™ adaptor set for RadioEar

Connects insert phones with Infant EarCups™. NB: Do not contain filters.  
Part no. Unit Price  
8500450 2 pcs £18.60

ABR

Easy ABR newborn hearing screening with the EARTurtle™

The EARTurtle™ disposables ear couplers is a new, soft, and easy-to-mount disposable ear coupler for ABR newborn hearing screening.  
You can use it with the RadioEar IP30 adapter and similar transducer couplers.

- Safe to apply on sensitive skin
- (complies with ISO 10993-1:2018)
- Consistent and reliable performance
- Durable and user-friendly design



Earturtle™ adapter for RadioEar IP30

Part no.	Unit	Price
8515872	1 pc	£8.03



Earturtle™ earcouplers

Part no.	Description	Unit	Price
8529735	w/tab electrodes & alcohol pads	36 sets	£282.40
8530442	w/snap electrodes & alcohol pads	20 sets	£156.42
8522749	w/o electrodes and alcohol pads	20 sets	£148.81

ABR LATEX



Sanibel™ Infant EarCup™ adaptor kit

Includes adaptors, 1 blue and 1 red tube.  
Part no. Unit Price  
8100577 Set/2 pcs £29.59

ABR





# Use Sanibel™ probe tubes for any equipment



The “REM effect” - better fitting & fewer follow-ups!



We often hear from the industry that only the manufacturers’ specific probe tubes should be used with their specific equipment. However, a recent comparative study by Sanibel Supply® disproves this myth.

In this study, our in-house audiologists and engineers compared the effects of different probe tubes on an Affinity2.0 system. They compared AudioScan tubes (red), Otometrics tubes (blue), Sanibel™ tubes (orange) and Sanibel™ tubes (grey).

This involved running each tube via the ‘tube calibration’ process and then comparing the following REUG curve when the tube had been left in the calibration position.

The investigation showed no significant difference in the REUG curve and all probe tubes achieved the same result when the calibration process had been run correctly.




Sanibel™ REM probe tubes, type A

Universal fit, length 76mm, beveled 45° angle tip and black marker ring.

Part no.	Unit	Price
8104116	Tray/50 pcs	£49.46

HIF



Sanibel™ probe tubes, type B

Fit all brands, **Interacoustics**, **Natus** Aurical, MedRx and AudioScan

Part no.	Unit	Price
8500360	Bag/40 pcs	£29.59

HIF



Silicone tubes for SPL60™

For use with Affinity™ REM probe.

Part no.	Unit	Price
8500700	Bag/50 pcs	£8.37

HIF



Sanibel probe tubes, type A, black, universal

Part no.	Unit	Price
8533301	packet/50 pcs	£31.28



Sanibel probe tubes, type C, blue, universal

Part no.	Unit	Price
8533303	packet/50 pcs	£31.28



Sanibel probe tubes, type D, blue, universal


Part no.	Unit	Price
8533304	packet/50 pcs	£31.28



Bionix Safe Ear Currettes

Safe Ear Currettes allows safe and gently removal of different types of cerumen. Single use ear currettes. Ensure your patients' safety by using these affordable ear currettes. Available as variety pack or single size items. Safe Ear Currettes allow healthcare provides to safely and gently remove cerumen with styles to match different types of wax and personal preference.

Part no.	Description	Unit	Price
8515585	Green MicroLoop®	box/50 pcs	£74.32
8515586	Orange ControLoop®	box/50 pcs	£74.32
8515587	Red AngleLoop®	box/50 pcs	£74.32
8519359	Purple VersaLoop®	box/50 pcs	£74.32
8519360	Blue InfantScoop®	box/50 pcs	£74.32
8519372	Yellow CeraSpoon®	box/50 pcs	£74.32
8519361	White FlexLoop®	box/50 pcs	£74.32
8519362	Variety Kit, (10 each style & 15 white)		£83.62



Bionix Lighted Ear Currettes

Box of 50 pcs. Use the Bionix lighted ear currettes for safe removal of earwax. The Bionix currettes comes in a variety of different shapes for you to option for your optimal choice in a situation. Bionix Lighted Ear Currettes combine illumination and magnification to help healthcare providers take less time to see and safely remove cerumen from the ear canal.

All versions including 1 Light Source and 1 magnifying lens.

Part no.	Description	Price
8519364	Lighted AngleLoop®, 4mm diameter	£93.01
8519365	Lighted MicroLoop®, 3mm diameter	£93.01
8519366	Lighted WaveCurette™, 2.5mm wide tip	£93.01
8519367	Lighted FlexLoop®, 4mm wide tip	£93.01
8519368	Lighted CeraSpoon®, 4mm wide tip	£93.01
8519369	Lighted InfantScoop®, 2mm wide tip	£93.01
8519373	Lighted VersaLoop®, 3mm wide tip	£93.01
8519370	Lighted CeraPik™, 2.5mm concave tip. Great for hard, impacted cerumen and foreign bodies.	£93.01



Lighted Ear Curette Variety Pack

Part no.	Description	Price
8519363	with 1 Light source and magnifying lens includes all tip styles except CeraPik™	£111.52



Lighttec

Bionix Lic illumination and magnification to help healthcare providers take less time to see and safely remove ear canal and nasal passage obstructions.

Part no.	Description	Price
8519392	10 foreceps, 1 light and 1 magnifying lens	£120.82



OtoClear Irrigation Tips

Bionix OtoClear Ear Irrigation helps provide safe and comfortable ear irrigation without the backsplash.

Part no.	Unit	Price
8519375	Box/40 pcs	£55.72



OtoClear Ear Basin

Part no.	Unit	Price
8519390	Box/3 pcs	£27.90



OtoClear Spray wash Kit

OtoClear Ear Irrigation system directs water to the ear canal wall gently with three divergent streams, instead of directly at the Tympanic Membrane

Part no.	Description	Price
8519371	box 20 irrigator tips, 1 470ml spray wash bottle, 1x ear basin	£92.92



O Scope

Part no.	Unit	Price
1000437	1	£1,375.00



Opticlar Vm5 Wireless Headlight With 2.5 x (340mm) Loupes

Part no.	Unit	Price
8535325	1	£1,375.00



1 Litre Gelsax Bacticlear Standard Liner = 1 Box of 25 Liners

Part no.	Unit	Price
90015824	1 box	£60.00



Suction Tubing 6mm x 3m

Part no.	Unit	Price
90015851	Pack/50 pcs	£95.00



Earol Olive Oil Spray 10ml

Part no.	Unit	Price
90015852	Pack/24 pcs	£140.00



Zoellner Suction Tube 30deg

Part no.	Unit	Price
90015853	Pack/50 pcs	£55.00



Zoellner 30 Deg W/Hole

Part no.	Unit	Price
90015854	Pack/50 pcs	£58.00



Zoellner 35 Deg Flat Handle

Part no.	Unit	Price
90015855	Pack/50 pcs	£55.00



Zoellner Fine End 18 SWG

Part no.	Unit	Price
90015856	Pack/50 pcs	£40.00





Jobson Horne Probes

Part no.	Unit	Price
90015857	Pack/100 pcs	£48.00



Wax Hooks

Part no.	Unit	Price
90015858	Pack/100 pcs	£50.00



Ear Specula Paediatric  
2.75mm

Part no.	Unit	Price
90015859	Pack/850 pcs	£50.00



Tumarkin Specula Size 6

Part no.	Unit	Price
90015866	Pack/100 pcs	£75.00



Tumarkin Specula Size 7

Part no.	Unit	Price
90015867	Pack/100 pcs	£75.00



Gruber Speculum No. 1

Part no.	Unit	Price
90015868	Pack/100 pcs	£75.00



Ear Specula Adult 4.25mm

Part no.	Unit	Price
90015860	Pack/850 pcs	£50.00



Tumarkin Specula Size 1

Part no.	Unit	Price
90015861	Pack/100 pcs	£75.00



Tumarkin Specula Size 2

Part no.	Unit	Price
90015862	Pack/100 pcs	£75.00



Gruber Speculum No. 2

Part no.	Unit	Price
90015869	Pack/100 pcs	£75.00



Gruber Speculum No. 3

Part no.	Unit	Price
90015870	Pack/100 pcs	£75.00



Gruber Speculum No. 4

Part no.	Unit	Price
90015871	Pack/100 pcs	£75.00



Tumarkin Specula Size 3

Part no.	Unit	Price
90015863	Pack/100 pcs	£75.00



Tumarkin Specula Size 4

Part no.	Unit	Price
90015864	Pack/100 pcs	£75.00



Tumarkin Specula Size 5

Part no.	Unit	Price
90015865	Pack/100 pcs	£75.00



Shea Aural Specula 4mm

Part no.	Unit	Price
90015872	Pack/100 pcs	£72.00



Shea Aural Specula 5mm

Part no.	Unit	Price
90015873	Pack/100 pcs	£72.00



Shea Aural Specula 6mm

Part no.	Unit	Price
90015874	Pack/100 pcs	£72.00



Croc Forceps  
Micro 3mm Sterile

Part no.	Unit	Price
90015875	1 pc	£8.50



Henckel Forceps  
80mm Sterile

Part no.	Unit	Price
90015876	1 pc	£8.50



Micro Scissor  
5mm Jaw Sterile

Part no.	Unit	Price
90015877	1 pc	£15.00



St Barts Blunt Hook  
Sterile

Part no.	Unit	Price
90015878	1 pc	£4.50



Formby Cerumen  
Scoop & Hook Sterile

Part no.	Unit	Price
90015879	1 pc	£4.50



Tilley Aural Forceps  
140mm Sterile

Part no.	Unit	Price
90015880	1 pc	£5.00



Stapes Hook  
Straight Sterile

Part no.	Unit	Price
90015884	1 pc	£6.00



Croc Forceps  
Std 5mm Sterile

Part no.	Unit	Price
90015885	1 pc	£8.50



Biopsy Forceps 2mm x 1mm,  
85mm Length Sterile

Part no.	Unit	Price
90015886	1 pc	£14.00



Biopsy Forceps 3mm x 2mm,  
85mm Length Sterile

Part no.	Unit	Price
90015887	1 pc	£14.00



Jobson Horne Probe  
178mm

Part no.	Unit	Price
90015888	1 pc	£2.00



Croc Forceps 80mm Shaft,  
5mm, Fine Jaw

Part no.	Unit	Price
90015889	1 pc	£8.50



Disp Ointment Syringe

Part no.	Unit	Price
90015881	Pack/50 pcs	£60.00



Stapes Hook Curved Sterile

Part no.	Unit	Price
90015882	1 pack	£6.00



Stapes Hook Angled Sterile

Part no.	Unit	Price
90015883	1 pc	£6.00



EarWay Pro -  
Wax Removal Helix

Very easy and safe way to remove ear wax. The helix grabs around the ear wax lump and pulls it out in a comfortable way

Part no.	Unit	Price
8518025	Pro M, 5,4mm, purple 25 units	£70.18
8518024	Pro S, 4,2mm, blue 25 units	£70.18
8521854	Pro L, 6,5mm, green 25 units	£70.18



PVC 30m Yellow  
Silicone tubing

Part no.	Unit	Price
90015726	1	£40.00



1 Litre Cannister

Part no.	Unit	Price
90015727	1	£20.00






Box of 25 1 litre liners

Part no.	Unit	Price
90015728	25	£40.00



1 Tapered Connector

Part no.	Unit	Price
90015729	1	£6.00



25m roll of 12mm ID silicon Tubing

Part no.	Unit	Price
90015713	1	£145.00



SAM 2 connection Tube with 2 x Elbow Connectors

Part no.	Unit	Price
90015718	3 pcs	£16.70



Vacsax 2 litre Cannister

Part no.	Unit	Price
90015719	1	£40.00



Vacsax 2 litre liners

Part no.	Unit	Price
90015720	Box of 30	£116.00



Box of 24 disposable Filters

Part no.	Unit	Price
90015710	24	£65.00



Box of 24 Hydrophobic Filters

Part no.	Unit	Price
90015711	24	£120.00



Pack of 10 Vacuum elbow Connectors

Part no.	Description	Price
90015712	Pack of 10	£35.00



Vacsax tapered connectors

Part no.	Unit	Price
90015721	Pack of 10	£65.00



Vacsax connection Tube with elbow connector and Vacsax tapered connector

Part no.	Unit	Price
90015722	3 pcs	£27.00



Replacement 2 litre SAM 2 autoclavable Jar

Part no.	Unit	Price
90015714	1 pc	£92.00



2m Silicon Patient Tube with fitted ends

Part no.	Unit	Price
90015715	1 pc	£10.00



SAM MS Unit for Microsuction

The SAM MS aspirator is a top quality surgical suction unit designed to run quietly for long periods of time.

At it's heart is a state of art oil-free diaphragm pump which is combined with multiple safety features including internal filters and overflow control valves.

- Vacuum 0 to 0.85 bar / 0 to 650mmHg
- Flow rate 30 l/min
- Collection Jars 1 x 2 litre
- Weight 10kg
- Noise Level 50dB


Part no.	Unit	Price
1000425	1	£750.00



Impression Pad, foam  
Standard Pink

Diameter approx 1.1cm, length 1.2cm


Part no.	Unit	Price
8515559	bag/100 pcs	£18.18



Foam Impression Pads,  
conical, Light Pink

Diameter approx 1.25cm


Part no.	Unit	Price
8515560	bag/50 pcs	£20.46



Impression Pad,  
foam, Large blue

Diameter approx 1.25cm

Part no.	Unit	Price
8515423	bag/50 pcs	£14.37



Impression Pads

Diameter approx. 1.1 cm, length 0,7 cm

Part no.	Description	Unit	Price
8522182	Foam Impression Pads CIC, yellow	100 pcs	£18.18



Impression Pads

Diameter approx. 1.1 cm, length 0,7 cm

Part no.	Description	Unit	Price
8527869	Foam Impression Pads Slim, extra small, white	50 pcs	£9.30



Mixing tray with  
dosing scale

Part no.	Unit	Price
8522105	1 pc	£13.95



Impression Pad, foam,  
Small Grey 9mm

Diameter approx. 0.95cm

Part no.	Description	Unit	Price
8515429	Grey,	bag/50 pcs	£14.37



Impression Pad, foam,  
Medium White 11mm

Diameter approx 1.1cm

Part no.	Description	Unit	Price
8515430	White	bag/50 pcs	£14.37



Cotton impression pad Small  
with blue safety thread

Diameter pprox 0.65cm

Part no.	Description	Unit	Price
8515571	Blue	bag/50 pcs	£16.06



Large Mixing Spatula


Part no.	Unit	Price
8522104	1 pc	£27.90



Impression Syringe Grey

Easy to fill, inject and clean. Provides a stable and easy injection. 4mm.

Part no.	Unit	Price
8515445	1 pc	£18.60



Mixing pad

Part no.	Unit	Price
8522106	1 pc	£7.44



Cotton Impression  
Pads, medium with  
white safety thread

Diameter approx 0.75cm

Part no.	Description	Unit	Price
8515572	White	bag/50 pcs	£16.06



Cotton Impression Pads,  
large with red safety thread

Part no.	Description	Unit	Price
8515573	Red	bag/50 pcs	£16.06



Cotton Impression Pads

Part no.	Description	Unit	Price
8527871	Cotton Impression Pads, extra small with white safety thread. Diameter approx 0.85cm	50 pcs	£11.41
8527870	Cotton Impression Pads, mini with white safety thread	50 pcs	£11.41



Otoform® Ak

Non-sticky, very pliable impression silicone. Final hardness approx. 35 Shore A

Note: requires an impression gun and mixing tip.

Part no.	Unit	Price
8530889	50ml double cartridges, 1 pc	£65.10



Otoform® Ak, Tubs (A+B),  
inc. dosing spoons

8528835 1box, (2x544ml) £60.45

Part no.	Unit	Price
8528835	1 box, (2x544ml)	£60.45



Fotoplast Lacquer 3

suitable for all hard earmolds. Curing time: 60seconds. Tissue-friendly surface sealing.

Part no.	Unit	Price
8530893	50ml	£32.13





### Manual Impression Gun

Item 8522434 to be used with 8522433.

Part no.	Unit	Price
8522434	Injection for 50ml double cartridges 1 pc	£73.56



### Fotoplast Gel

Light curing repair resin.

Part no.	Description	Price
8530895	clear transparent, 20ml	£32.13
8530896	opaque H (yellowish), 20ml	£32.13



### Otoform® A flex

Soft, flexible, easy processing. For CIC impressions, narrow auditory canals or people sensitive to pain. Final hardness approx. 30Shore A. Note: requires impression gun and mixing tip.

Part no.	Unit	Price
8522433	50ml double cartridges, 1 pc	£65.10



### Mixing Tips for 50ml double cartridge

For use with 8522433

Part no.	Unit	Price
8522432	5,4mm, Bag/100 pcs	£60.45
8522183	6,5mm, Bag/100 pcs	£60.45



### Impression material A-ZOFT

Part no.	Unit	Price
8515424	2 x 800g jars	£76.94
8515426	40 foil bags (sets 2 x 8g)	£69.75



### Bite Blocks

Disposable styrofoam spacers hold a patient's mouth open while taking open jaw impressions.

Part no.	Unit	Price
8527849	bag/100 pcs	£51.15
8527850	bag/25 pcs	£23.25



### Cerumen Removal Tool

Part no.	Description	Unit	Price
8527907	with Brush & Magnet	1 pc	£2.79
8527908	with Brush	1 pc	£1.86
8527909		1 pc	£1.35



### LED Earlite Tips

To be used with LED Earlite 8527868

Part no.	Unit	Price
8527862	LED Earlight Straight Tip Bag/4 pcs	£9.30
8527863	LED Earlight Pros Tip Bag/4 pcs	£12.09
8527864	LED Earlight Scoop Tip Bag/4 pcs	£8.37
8527865	LED Earlight Steel loop Tip Bag/2 pcs (small & large)	£21.39



### LED Earlite

The 18 Lumens light is concentrated in a small tip that fits most ears. Includes 2 AAA batteries

Part no.	Unit	Price
8527868	1 pc	£13.95



### Battery Extractor Tool

Part no.	Unit	Price
8527911	1 pc	£7.44



### Vent Cleaning Tool

Part no.	Unit	Price
8527910	1 pc	£1.86



### Mini Cleaning Brush

Clean hearing instruments in an easy way.

Part no.	Unit	Price
8519100	1 pc	£1.10



### Earlight with Scaled Top

Part no.	Description	Price
8519357	Luxamed Ear light black, acrylic scaled top	£23.25
8519358	Luxamed Acrylic Top with Scale	£1.69



### Stethoclip

Hearing Instrument listening stethoscope for reasuring the fitting process

Part no.	Description	Unit	Price
8515600	Stethoclip, black	1 unit	£18.60
8520346	Ear funnels	1 unit	£2.79
8520347	ITE Adaptor	1 unit	£5.16

## Take care of your hearing aids

Take good care of your hearing aid. Your Hearing instrument is a high tech device needing regular care in order to provide you the best performance every day.



### Cleansing tablets

Part no.	Unit	Price
8520148	4 pk	£3.30



### Dehydrating capsules

Part no.	Unit	Price
8520149	4 pk	£3.30



### Drying cup

Part no.	Unit	Price
8520150	1 pc	£1.44



### Cleaning cup

Part no.	Unit	Price
8520151	1 pc	£1.44



### Otofloss

Part no.	Unit	Price
8520152	1 pc	£10.99



### Perfect Dry LUX box for all types of hearing aids

Part no.	Unit	Price
8520154	1 pc	£55.72



### Starter kit

Part no.	Unit	Price
8520155	kit	£10.99



### Drying spray

Part no.	Unit	Price
8520153	1 pc	£6.51



### Earmold Air Blowers

Part no.	Unit	Price
8515562	Box/12 pcs	£3.30



### Cleaning spray

Part no.	Unit	Price
8520147	1 pc	£3.64



### Cleaning wipes

Part no.	Unit	Price
8520146	30 pk	£4.23



### Protect your ears

Part no.	Description	Price
8520052	Natural lubricant	£6.43



### Protect your ears

Part no.	Description	Price
8520053	Dermoprotective Gel	£4.57



### Protect your ears

Part no.	Description	Price
8520063	Earspray	£9.22



### Sanohra Earplugs Fly

Sanohra Fly ear plugs. Reduces ear pressure experiences when flying.

Part no.	Description	Unit	Price
8519351	EarPlugs Fly Kids	1 set	£8.88
8519352	EarPlugs Fly Adult	1 set	£8.88



### Sanohra Earplugs Swim

Sanohra Swim earplugs. Enables people to swim without getting water in the ear, reducing ear infections.

Part no.	Description	Unit	Price
8519353	EarPlugs Swim Kids	1 set	£7.02
8520143	EarPlugs Swim Adult	1 set	£7.02

### SANOHRA Fly ear plugs

Beyond "noise protection" it can be important to protect the ear drum (tympanic membrane) from changes in external pressure that are too rapid, e.g., during take off and landing phases in airplanes. For airplane passengers experiencing a problem with the pressure equalisation, the use of a special ear protection product means freedom from an unpleasant feeling or piercing pain during flight.

### SANOHRA Swim ear plugs

In case of injury to the ear drum or insertion of a vent tube, it is necessary to prevent water from reaching the ear drum through the outer auditory canal. Swimming, bathing or showering might cause inflammation of the middle ear. Staying in chlorinated water for an extended period of time can cause the skin in the outer auditory canal to swell and cause inflammation (otitis externa).



### Pre-bent tubing

Quilled on one end for easy insertion into the earpiece and double bent to lie close to the head

Part no.	Description	Unit	Price
8527852	Standard #12	25 pcs	£9.30
8527872	Medium #13	25 pcs	£9.30
8527853	Thick #13	25 pcs	£9.30
8527854	Extra Thick #13	25 pcs	£9.30



### Pre-bent Disappear Tubing

With a long quill for easy insertion into earmold. The tubing is perfectly formed, free from kinking, in flesh colour. Please note, tubing color is not guaranteed. Colour of tubes may vary from pictures

Part no.	Description	Unit	Price
8527855	#12S	10 pcs	£14.88
8527856	#13M	25 pcs	£20.46
8527857	#13T	25 pcs	£20.46



### Tubing Retention System (TRS)

Designed for all soft earpiece materials. Each TRS tube uses an acrylic ring to hold it in position in the earpiece. TRS tubes are easily replaced and will not tear silicone materials.

Part no.	Description	Unit	Price
8527927	Tubing Retention System#13T.	25 pcs	£27.90
8527928	Tubing Retention System#13ST	25 pcs	£27.90



### Superglue, viscous

Instant adhesive - low viscosity

Part no.	Unit	Price
8530962	bag/10 pcs	£17.76



### Arc Tubing connection

with retentions

Part no.	Unit	Price
8530938	bag/10 pcs	£19.45



### Tube Lock without tubing

Brass Tube-Locks Tube-Locks are used for attaching tubing to flexible type earmolds without requiring any cement, yet allowing tubing replacement as needed.

Part no.	Description	Unit	Price
8527858	#13 Thick	1 pc	£1.10
8527859	#13 Super Thick	1 pc	£1.10



### Pre-bent Tubing with Tube Lock

Tubing system that uses a brass retention ring for soft earpieces.

Part no.	Description	Unit	Price
8527860	#13 Thick	10 pcs	£18.60
8527861	#13 Super Thick	10 pcs	£18.60



### Pre-bent Stay Dry Tubing

Specially formulated tube helps prevent moisture buildup caused by humidity or perspiration

Part no.	Description	Unit	Price
8527873	Pre-Bent Tubes Stay Dry #13M	10 pcs	£17.67
8527874	Pre-Bent Tubes Stay Dry #13T	10 pcs	£17.67



### Tubing

Ø 2.0 x 3.0mm

Part no.	Unit	Price
8530939	1 pc/10 meters	£21.14



### Tubing Floss

For cleaning tubing

Part no.	Unit	Price
8530957	OtoVita Tubing Floss micro tubing, 30 pcs	£12.68
8530958	OtoVita Tubing Floss tubing, 30 pcs	£12.68



### Horn Tube

Special tube designed to increase gain over conventional tubing. Ideal when using an analog hearing instrument and fitting a high frequency hearing loss. Also designed to allow for a friction fit into earmold without the need to damage soundbore when removing it. Diappear Colour

Part no.	Description	Price
8515663	TRS Tube, 4mm, Dry-Tube	£45.57
8515657	TRS Tube, 3mm, Dry-Tube	£45.57



### Horn Tube

Uses the horn effect to increase response in the high-frequency range

Part no.	Description	Unit	Price
8527934	Horn Tub 3mm Medium	6 pcs	£27.90
8527930	HornTube 3mm Thick	6 pcs	£27.90
8527931	Horn Tube 4mm Medium	6 pcs	£27.90



### Horn TRS Tubes

Horn TRS tubes use acrylic rings to hold tubes in position. Designed for soft earpiece materials.

Part no.	Description	Unit	Price
8527932	Horn TRS Tube w/acrylic ring3mm	6 pcs	£32.55
8527933	Horn TRS Tube w/acrylic ring4mm	6 pcs	£32.55







### Surgical Blades & Handles

Suitable for cutting HA tubing.

Part no.	Description	Price
8515427	Surgical blade handle No.3	£10.23
8515428	Surgical blade No.11, 100 pcs	£27.65



### Loctite 406 Instant adhesive

For fast bonding of plastics and elastomeric materials. This product offers exceptionally low viscosity.

Part no.	Unit	Price
8005067	1 unit	£18.60



### Probe Set, Stainless steel

For probing, material removal.

Part no.	Unit	Price
8527929	set of 4 pcs	£27.90



### Scissors, 3.5"

Sharp scissors suitable for general use in the clinic.

Part no.	Unit	Price
8527905	1 pc	£6.51



### Tubing Expander Plier

Stainless steel pliers with two-way lever spring allows you to quickly and easily stretch acoustic tubing before fitting it over earhooks, elbows, etc.

Part no.	Unit	Price
8527875	Curved, 1 pc	£41.85
8527876	Straight, 1 pc	£41.85



### Tube Puller, Monofilament

Soft and flexible enough to easily slip through the sound bore of an acrylic earpiece, allowing quick and simple tube replacement.

Part no.	Unit	Price
8527877	1 pc	£12.09



### Tube Puller, Wire Loop

Rigid enough to slip easily through the sound bore of a vinyl or silicone earpiece, allowing quick and simple tube replacement.

Part no.	Unit	Price
8527878	1 pc	£13.02



### End Nipper Pliers

A perfect tool for trimming tubing from the end of a canal after the earmold has been retubed.

Part no.	Unit	Price
8527900	1 pc	£12.09



### Needle Nose Pliers 5"

Needle nose pliers for a multitude of uses in the hearing healthcare clinic.

Part no.	Unit	Price
8527898	1 pc	£12.09



### Sprue Cutter, 4.5"

A unique design easily trims tubing flush to the earpiece.

Part no.	Unit	Price
8527935	1 pc	£12.09



### Superglue, viscous

Instant adhesive - low viscosity

Part no.	Unit	Price
8530937	1 pc/10ml	£10.15



### SK 87 Silicone glue

Part no.	Unit	Price
8532585	1 pc/18ml (20g)	£11.84



### Alligator Forceps

Standard 5.5" stainless steel alligator forceps with serrated jaw. Not intended for invasive medical procedures

Part no.	Unit	Price
8527901	1 pc	£51.15



### Tweezers

For a variety of clinical uses. Good for working on hearing instruments and other small parts.

Part no.	Unit	Price
8527902	Curved Point Precision Tweezers 6", 1 pc	£11.16
8527903	Fine Tipped 4.75" Tweezers, 1 pc	£13.95



### Precision Screwdriver set

For making adjustments to hearing instruments or other precision work.

Part no.	Unit	Price
8527896	set of 6 pcs, Includes four flat-tip screwdrivers and two Phillips tips.	£7.44
8527897	set of 11 pcs, Six flat-tip screwdrivers, three Phillips tips, one needle point tip, one magnet tip.	£11.16

# Surface electrodes tailored to your need

We believe that it is important to choose the correct surface electrode for each specific recording and patient type.

### Wetgel or hydrogel electrodes?

A wetgel electrode consists of water, thickener, ionic salts, surfactants and bactericides or fungicides. Wetgel migrates easily and rapidly onto the skin surface, reducing impedance levels, which in turn will generally result in lower noise in the recording. A hydrogel surface electrode is less adhesive when compared to wetgel, making it particularly suitable for newborns.

### Chloride salt concentration

The higher the chloride salt content, the more conductive the electrode. Higher salt content (10% and above), pre-gelled surface electrodes are useful for making fast, high quality measurements of biopotentials. High conductivity electrodes generally have reduced artifact, due to the lower impedance generated between electrode and skin surface. As the chloride salt content drops, there is a corresponding reduction in electrode conductivity.

If the chloride content is below 10%, the electrode can then be employed for longer recording times, with a reduced risk of skin irritation. As a rule of thumb, wetgel has a chloride content of 10% and hydrogel 4%.

It is important to use electrodes with similar chloride salt content when testing.

The impedance of the electrode/skin junction is highly dependent on the electrolyte type and the chloride salt concentration. For example, a hydrogel electrode with 4% chloride concentration will have approximately 10x higher impedance than a wet liquid gel electrode with 10% chloride concentration, after first application to the skin.

### Adhesion

An advantage of a conductive adhesive is that the active surface can be as large as the electrode itself because no additional adhesive strip is needed. The Sanibel disposable electrodes come with a range of adhesive qualities. All electrodes are designed to adhere well to skin surfaces, but the adhesion characteristics can be grouped according to the application:

- Strong adhesive
- Moderate adhesive

Strong adhesion is not recommended for newborns.

### Skin preparation

For optimum electrode to skin conductivity, the skin should be lightly brushed with a gentle abrasive gel/wipe. An alcohol wipe is not recommended to improve conductivity as this will only dry out the skin surface.



### Sanibel™ snap electrodes

Pre-gelled with a very gentle hydrogel, perfect for ABR testing on infants. Size 2.2 x 3 cm. Item no. 8107137 is packed in 20 smaller pouches of 3 electrodes.

Part no.	Unit	Price
8107137	Box/60 pcs (20x3 pcs)	£35.09

ABR



### Sanibel™ tab electrodes

Pre-gelled cloth electrode with a very gentle hydrogel, perfect for ABR testing on infants. Size 2.2 x 3 cm. Packed in 20 smaller pouches of 3 electrodes.

Part no.	Unit	Price
8107410	Box/60 pcs (20x3 pcs)	£27.82

ABR



### Vermed Hydrogel Snap Electrode

Large soft electrode size suitable for newborn screening and ABR measurements. Cloth material makes it convenient to wear. Size: 22.2. Intended use areas: Newborn, ABR.

Part no.	Unit	Price
8516907	Bag/60 pcs (3x20 pcs)	£16.49

ABR



### E-A-RLINK™ gold foil tip trodes

For ECoG. Use with cable 8011343.

Part no.	Description	Unit	Price
8500370	Adult 13mm	20 pcs	£44.64
8500440	Pediatric 10mm	20 pcs	£44.64

ABR



### Sanibel™ TM electrode

For ECoG. Use with cable 8500650.

Part no.	Unit	Price
8011489	Box/10 pcs	£186.01
8049550	Single pcs	£22.32

ABR




### Vermed Snap Solid Gel Z35

Solid Gel type snap electrode suitable for newborn screening and Diagnostic follow-up. Latex free. Intended use areas: Newborn, Diagnostic follow-up all ages.

Part no.	Unit	Price
8516252	Bag/50 pcs	£6.34

ABR



### Ambu® Ground electrode 71410-M/1

Part no.	Unit	Price
8500680	1 pc	£1.86

ABR



### Vermed Tab Electrode

Rectangular tab electrode for newborn screening. 12,7x22,9mm. Conductive Adhesive Hydrogel. 10x10 pcs. per pouch. Intended use areas: Newborn, O-Vemp.

Part no.	Unit	Price
8516225	Bag/100 pcs	£8.79

ABR



### Vermed Wetgel Snap Electrode

Pre-gelled hydrogel electrode with Ag/AgCl sensor. Suitable for newborn screening or diagnostic follow up. Intended use areas: Newborn, Diagnostic follow-up, ABR, C-Vemp.

Part no.	Unit	Price
8516221	Box/75 pcs	£27.90

ABR





**Vermed Hydrogel Snap Electrode Square**

Smaller than the rectangular version, but otherwise similar. Good soft electrode size making it very suitable for newborn screening and ABR measurements. Cloth material makes it convenient to wear. Intended use areas: Newborn, ABR.

Part no.	Unit	Price
8102022	60 pcs (3x20 pcs)	£16.49

ABR



**Vermed Tab Electrode**

Rectangular tab electrode for newborn screening. 22.9x34mm. Conductive Adhesive Hydrogel. 10x10 pcs. per pouch. Intended use areas: Newborn, O-Vemp.

Part no.	Unit	Price
8516227	Baq/100 pcs	£6.09

ABR



**Vermed Economy line Snap Wet Gel Electrode**

Part no.	Unit	Price
8516251	Baq/50 pcs	£7.61

ABR



**Ambu® 0415M tab electrode**

Medium size tab electrode with solid gel for 12-lead ECG. Presented conveniently with 10 pieces to a card. Fish tail tab design reduces electrode lift-offs. Suitable for both adults and children. Intended use areas: ABR Diagnostics of all kinds with Tab electrodes.

Part no.	Unit	Price
8516138	Baq/100 pcs (10x10 pcs)	£10.23

ABR



**Ambu® White Sensor 40554**

Part no.	Unit	Price
8520431	Baq/3 pcs	£3.30

ABR



**Ambu® Blue Sensor™ N-50-K**

Foam ECG electrode with 10 cm wire and 4mm connector, for use with universal electrode cable. Wetgel.

Part no.	Unit	Price
8521842	Baq/25 pcs	£25.37

ABR



**Ambu® 2300 Blue Sensor Tab electrode**

The medium size electrode features a skin-friendly solid gel providing optimal signal quality recordings during 12-lead diagnostic resting ECG application. Suitable for ABR diagnostics of all kinds. Intended use areas: ABR Diagnostics of all kinds.

Part no.	Unit	Price
8516139	Baq/100 pcs	£13.11

ABR



**Ambu® Neuroline 720**

Pre-gelled and self adhesive, ideal for use in sleep registration where the strong, yet skin-friendly, adhesion ensures that the electrode stays in place for the entire examination, even if the patient sweats. Intended use areas: ABR Diagnostics of all kinds.

Part no.	Unit	Price
8500060	Baq/25 pcs	£12.09

ABR



**Ambu® Neuroline 720 with 100cm cable & k-connector**

Pre-gelled and self adhesive, ideal for use in sleep registration where the strong, yet skin-friendly adhesion ensures that the electrode stays in place for the entire examination, even if the patient sweats. Intended use areas: ABR Diagnostics of all kinds.

Part no.	Unit	Price
8503244	Baq/10 pcs	£16.74

ABR



**Ambu® Blue Sensor™ N-00-S**

Foam ECG electrode for use with button cables. Wetgel.

Part no.	Unit	Price
8013029	Baq/25 pcs	£13.53

ABR



**Ambu® Blue Sensor™ P-00-S**

High quality Ag/AgCl sensor - ensures excellent trace quality.

Part no.	Unit	Price
8504782	Baq/50 pcs	£14.88

ABR



**Ambu® Blue Sensor™ N-10-A**

Foam ECG electrode with 10 cm wire and 4mm connector, for use with universal electrode cable. Wetgel.

Part no.	Unit	Price
8500670	Baq/25 pcs	£21.39

ABR

Chosing the correct electrode for optimal performance

Interacoustics Supplies Hub offers various kinds of surface electrodes for ECG as well as EEG and ABR. Surface electrodes vary in size, adhesion, connectivity and frequency (*test category*).

Our surface electrodes are always tested and validated to ensure optimal performance for your testing!





Sanibel™ ECoChG starter kit

For Eclipse™ and Bio-logic® AEP system. Includes cable, gel and 2 pcs TM electrodes.

Part no.	Unit	Price
8013209	Set/4 pcs	£93.01

ABR



Pinch Clip electrode cables 100cm

Part no.	Unit	Price
8519048	3 pcs	£35.34



Sanibel™ ECoChG cable

For Eclipse™ and Bio-logic® AEP system. 1.5mm/.060" female connector. For use with Sanibel™ TM electrodes.

Part no.	Unit	Price
8500650	1 pc	£80.32

ABR



Tip trode cables

For use with TTE25 gold foil tip trodes. 1.5mm with touch proof connectors and plated clamp clips. 50 cm long.

Part no.	Unit	Price
8011343	Set/2 pcs	£58.17

ABR



Ambu electrode cable, 76cm/30"

Allows two of the preamplifier inputs to be connected, e.g. for linking mastoid electrodes for cortical testing.

Part no.	Unit	Price
8518508	1 pc	£14.88

ABR




Sanibel™ EP electrode jumper cable

Allows two of the preamplifier inputs to be connected, e.g. for linking mastoid electrodes for cortical testing.

Part no.	Unit	Price
8107872	1 pc	£14.88

ABR



Disc electrodes

For use with Tenzo™ electrode paste. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). 50 cm cable.

Part no.	Unit	Price
8011498	Set/3 pcs RWB	£69.75
8011345	Set/5 pcs RWBYK	£97.66

ABR




Universal electrode cable

Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). 50cm cable. Non-tangling.

Part no.	Unit	Price
8011353	Set/3 pcs RWB	£35.34
8011344	Set/5 pcs RWBYK	£58.59

ABR




Reusable golden-plated cup electrodes

Red (R), White (W), Blue (B), Yellow (Y), Black (K)

Part no.	Unit	Price
8107418	Set/3 pcs RWB	£79.05
8107145	Set/5 pcs RWBYK	£111.61

ABR



Alligator clip cable

For use with tab electrodes. 1.5mm/.060" female connector. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K).

Part no.	Unit	Price
8516581	Set/3 pcs RWB	£38.05
8513776	Set/3 pcs WYK	£38.05
8516584	Set/5 pcs RWBYK	£63.41

ABR



Pinch style cable

For use with snap electrodes. 1.5mm/.060" female connector. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K).

Part no.	Unit	Price
8500390	Set/3 pcs RWB	£35.34
8500400	Set/3 pcs WYK	£35.34
8500410	Set/5 pcs RWBYK	£58.59

ABR



Standard button cable

For use with snap electrodes. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). 50 cm cable. Non-tangling.

Part no.	Unit	Price
8532948	Set/3 pcs RWB	£38.05
8532949	Set/5 pcs RWBYK	£63.41

ABR



**SignaGel®, 250g**  
Hypoallergenic, bacteriostatic  
conductive saline electrode gel.

Part no.	Unit	Price
8102628	Tube/250g	£4.65

ABR



**Spectra 360®  
electrode gel, 250g**  
Salt- & chloride-free conductive  
electrode gel suitable for  
long-term applications.

Part no.	Unit	Price
8500740	Tube/250g	£5.58

ABR



**Ambu® 2121M Skin Prep Pads**  
Efficient, yet gentle skin abrasive.  
Chloride and acetone free. Pack of  
1x4 oz/114 g tube or 3x4 oz/114 g tube.

Part no.	Unit	Price
8523414	Box/500 pcs	£37.20



12 rolls

**3M™ Transpore™  
Surgical tape, 9.14m**  
Transparent, latex-free, water-  
resistant and porous for breathability.  
25mm x 9.14m.

Part no.	Unit	Price
8500750	Box/12 rolls	£26.97

ABR ~~LATEX~~



**Tenzo™ conductive  
EEG paste, 114g/228g**  
Adhesive gel for disc electrodes.  
Available in pack of 3x4oz/114g  
tubes, 3x8oz/228g jar.

Part no.	Unit	Price
8012924	Box/3 tubes	£26.04
8500330	Box/3 jars	£41.85

ABR



**Everi™ abrasive paste, 160g**  
Reduces the skin impedance  
during EEG, EP, ECG and Holter  
Stress test examination.

Part no.	Unit	Price
8012952	Bottle/160g	£7.95

ABR



12 rolls

**3M™ Micropore™  
Surgical tape, 9.14m**  
Latex-free and hypoallergenic. Holds  
well on damp skin. 25mm x 9.14m.

Part no.	Unit	Price
8500760	Box/12 rolls	£23.25

ABR ~~LATEX~~



**B Braun®  
alcohol pads**  
70% isopropyl alcohol.

Part no.	Unit	Price
8500050	Box/100 pcs	£7.02

ABR OAE



**Sonavelle® electrode  
gel, 250ml**  
Formaldehyde-free and glycerol-free.

Part no.	Unit	Price
8500380	Bottle/250ml	£7.44

ABR



**Gauze swabs**  
Non-sterile, 75 x 75cm and 8 layered.

Part no.	Unit	Price
8500040	Pack/100 pcs	£4.65

ABR OAE



**3M™ Red Dot™  
trace prep, 5m**  
One-step skin preparation.  
18mm x 5 m / 0.70" x 196".

Part no.	Unit	Price
8500030	Dispenser/5m	£41.85

ABR



10  
packs

**Q-tips**  
Cotton tipped applicators for  
patient preparation. Non-sterile.  
Each pack contains 100 pcs.

Part no.	Unit	Price
8500810	Box/10 packs	£27.90

ABR OAE



**NuPrep™ skin prep gel, 114g**  
Efficient, yet gentle skin abrasive. Chloride and acetone free.

Part no.	Unit	Price
8500320	Box/3 tubes	£30.27

ABR





500 pcs

### 3" Disposable ear cushion cover

Replaces disinfectants. Headphone Covers for classic audiometric headphones.

Part no.	Unit	Price
8107419	Box/500 pcs	£63.41

AUD LATEX



### Probe cleaning tool

For **Interacoustics** Titan™/ AT235™/AA222™/Sera™/MT10™/ **Amplivox** Otowaves/**MAICO** Mi XX + touchTym™ and **GSI** TymStar™ Pro.

Part no.	Unit	Price
8011351	1 pc	£2.37

ABR IMP OAE



### 3-in-1 floss

Used to clean probe tips.

Part no.	Unit	Price
8003565	Box/100 pcs	£9.30

ABR IMP OAE



100 pcs

### 4.5" Disposable ear cushion cover

Replaces disinfectants. Headphone Covers for noise enclosure audiometric headphones. For HDA200 /DD450/DD65v2/Peltor headsets.

Part no.	Unit	Price
8513404	Baq/100 pcs	£19.45

AUD LATEX



### Disinfectant detergent, 1000ml

Alcohol-free, ready to use solution for a quick and extensive disinfection of non-invasive medical devices.

Part no.	Unit	Price
8501706	Bottle, 1000ml	£27.90



### OtoVita® Professional

Disinfection Concentrate, 1l

Part no.	Unit	Price
8530888	1 pc	£41.43



### Keyboard protective foil

Place a protective foil over your keyboard to minimize contamination with grime, bacteria and germs, or liquid.

Part no.	Unit	Price
8500460	1 pc	£13.95



### Bridge & Implant floss

Used to clean probe tips.

Part no.	Unit	Price
8003566	Box/30 pcs	£3.72

ABR IMP OAE



### Disinfecting Wipes

Specially designed for surfaces, instruments, healthcare seats and medical devices.

Part no.	Unit	Price
8525161	Canister/90 wipes	£10.15



### Nitrile examination gloves

For use at medical sites and at research institutes etc. Do not interact with petroleum oils, grease, solvents. Do not contain latex.

Part no.	Description	Unit	Price
8515551	Small	100 pcs	£18.00
8515552	Medium	100 pcs	£18.00
8515553	Large	100 pcs	£18.00
8515554	X Large	100 pcs	£18.00

### CE certified Safety goggles

Part no.	Unit	Price
8522417	1 pc	£7.02

### Disposable head covers

Part no.	Unit	Price
8519350	Box/100 pcs	£21.39

### Disposable apron

Part no.	Unit	Price
8519348	Box/100 pcs	£29.76

### Disposable shoe covers

Part no.	Unit	Price
8519349	Box/100 pcs	£18.60

A clean clinic is also part of your professional image

Infection control is so much more than cleaning probe tips, protecting your personal computer's keyboard, and maintaining high standards within hygiene using e.g. disinfectant detergent for a quick and extensive disinfection of non-invasive medical devices.

When it comes to infection control, preventing germs, dirt and dust from spreading via other sources, is also of critical importance - both to the safety of your patients and yourself! Gloves, shoe covers, facemasks, and aprons are just some of the "accessories", which can help you succeed in this.

Maintain high, professional standards and make your clinic an ever-safe environment for your patients, making the wearing of gloves, facemasks and other body protection wear a natural part of the "patient interaction" process.

How to clean probe tips!

Watch our video at:  
[www.interacoustics.co.uk/servicehub](http://www.interacoustics.co.uk/servicehub)





USB cables A to B

Part no.	Description	Unit	Price
8011241	2m/6.56 feet	1 pc	£13.95
8011243	5m/16.40 feet	1 pc	£20.21

AUD



TDH39™/DD45™ replacement cable, symmetric

For TDH39™ and DD45™. 30° jack plug.

Part no.	Description	Unit	Price
8507203	2m	1 pc	£32.13
8507204	5m	1 pc	£42.28

AUD



HDA200 replacement cable, symmetric

For HDA200 and DD65 headphones. 2 m cable and 6.3mm/ 30° jack plug.

Part no.	Unit	Price
8004172	1 pc	£32.13

AUD



Replacement cables for audiometry headsets, Single Mono 2m Contra Earphone Cable

Part no.	Unit	Price
8507205	1 pc	£32.13

AUD

IMP



Replacement cable for Bone Conductor B71/B81 (with safety fixation point)

Part no.	Unit	Price
8515713	1 pc	£32.13

AUD



IP30 Contra Insert earphone cable

Cable length: 2 m. Plug type: 30°, 1x6.3mm gold-plated mono jack.

Part no.	Unit	Price
8107359	1 pc	£32.13

AUD



RadioEar Spring Steel headband (P3333)

Part no.	Unit	Price
8507213	1 pc	£29.59

AUD

HIF

ABR



B71 Bone conductor replacement cable

2 m cable and 6.3mm/ 30° jack plug.

Part no.	Unit	Price
8508473	1 pc	£32.13

AUD



Peltor® headset/earmuff replacement hygiene kit

For DD450, HDA200 and DD65.

Part no.	Unit	Price
8500270	1 set	£26.02

AUD



RadioEar Foldable audiometric headband

Black, foldable.

Part no.	Unit	Price
8010882	1 set	£66.37

AUD



RadioEar RE-7 Headband

Black.

Part no.	Unit	Price
8517379	1 set	£49.88

AUD



RadioEar DD45 Adult audiometric headset

Headband: HB3045 .Transducer: DD45. Cable length: 2 m. Plug type: 30°, 2x6.3mm gold-plated mono jack.

Part no.	Unit	Price
8105268	1 set	£420.00

AUD



RadioEar DD45 Audiometric headset

Headband: HB7. Transducer: DD45. Cable length: 2 m. Plug type: 30°, 2x6.3mm gold-plated mono jack.

Part no.	Unit	Price
8010954	1 set	£420.00

AUD



RadioEar DD45 Adult audiometric headset

Headband: HBA, foldable. Transducer: DD45. Cable length: 2 m. Plug type: 30°, 2x6.3mm gold-plated mono jack.

Part no.	Unit	Price
8106340	1 set	£420.00

AUD





### RadioEar DD45 transducer

Audiometric earphones engineered to match the performance and function of the most popular earphones available for audiological testing. 10 ohm receiver. Easy installation. Reliable. RoHS compatible.

Part no.	Unit	Price
8507268	1 set	£54.11

AUD



### Headband cushion

For Sennheiser HDA200, RadioEar DD450, DD65v2 and Peltor audiometric headbands.

Part no.	Unit	Price
8000770	1 pc	£15.90

AUD



### Sanibel™ TDH39/DD45 earphone cushion red & blue, model 51

For TDH39, 49, 50 & DD45 earphones.

Part no.	Unit	Price
8107611	1 Set	£21.56

AUD



### Ear cushion for DD45 & TDH39

Black, Model 51.

Part no.	Unit	Price
8507206	1 pc	£11.84

AUD



### RadioEar DD65 v2 Passive attenuation audiometric headset

Cable length: 2 m. Plug type: Straight, 2x6.3mm gold-plated mono jack.

Part no.	Unit	Price
8520330	1 set	£420.00

AUD



### Foam pad for contra headset

Compatible with e.g. TDH39c™ and DD45c™.

Part no.	Unit	Price
8507292	Bag/10 pcs	£11.58

AUD



### RadioEar replacement pads (3288-1)

For RadioEar Spring Steel headband.

Part no.	Unit	Price
8504777	Bag/10 pcs	£16.06

AUD HIF ABR



### Ear cushion HZP 09 pair

Ear cushion foam (1 pair) suitable for: PC 20, PC 25, PC 30, PC 31, PC 35, PC 36, PC 130, PC 131, PC 135, PC 136, PC 140, PC 141, PC 145, PC 146.

Part no.	Unit	Price
8104416	1 pair	£7.02

AUD



### Headband cushion

For RadioEar HBA foldable audiometric headband.

Part no.	Unit	Price
8010886	1 pc	£4.65

AUD



### RadioEar AMBAND bone conductor headband (P1188)

For adult, children and infants.

Part no.	Unit	Price
8507214	Bag/10 pcs	£40.58

AUD HIF ABR







24 pcs

### Disposable VNG goggle foam pads

For **Interacoustics®** VNG mask.

Part no.	Unit	Price
8011527	Box/24 pcs	£27.90

BAL



### Disposable specula for Difra Airstar

Part no.	Description	Unit	Price
8011457	2.5mm spec. (green)	100 pcs	£14.80
8011456	4mm spec (blue)	100 pcs	£14.80



### Paper Kidney tray

Part no.	Description	Price
8520139	10 pcs	£1.48



34 pcs

### Disposable ear specula

For AirFx™.

Part no.	Description	Unit	Price
8013013	2.75mm	34 pcs	£6.09
8013014	4.25mm	34 pcs	£6.09

BAL



50 pcs

### Water irrigator tips

For use with **Interacoustics®** Aqua Stim™.

Part no.	Unit	Price
8013018	Bag/50 pcs	£37.20

BAL



### Thermal paper for ERO-SCAN™/OtoRead™, 25m

Width: 57mm. Inner diameter: 12mm  
Outer diameter: 45mm

Part no.	Unit	Price
8500200	1 roll/25m	£3.72

OAE



### AP1300 thermal paper for Titan™, 10m

Width: 57mm. Inner diameter: 12mm  
Outer diameter: 35mm

Part no.	Unit	Price
8500230	1 roll/10m	£3.72

ABR

IMP

OAE



### Thermal paper for Sanibel™ MPT-II printer, 20m

For the Sanibel™ MPT-II printer which is supplied with **MAICO®** ERO-SCAN™, **Grason-Stadler®** Corti™ and **Interacoustics®** OtoRead™.  
Width: 57mm. Inner diameter: 9mm  
Outer diameter: 40mm

Part no.	Unit	Price
8029305	1 roll/20m	£3.72

IMP

OAE



### TPR10™ thermal paper for MTP10™, 25m

Width: 79mm. Inner diameter: 12mm  
Outer diameter: 45mm

Part no.	Unit	Price
8500220	1 roll/25m	£3.72

IMP



### Sanibel™ MPT-III thermal paper, 20m

For Sanibel™ MPT-III printer supplied with **Interacoustics®** AT235™.  
Width: 79mm. Inner diameter: 12mm  
Outer diameter: 45mm

Part no.	Unit	Price
8500610	1 roll/20m	£3.72

ABR

IMP

OAE



### TPR26™ thermal paper for Interacoustics®, 25m

Fits most popular audiometers and tympanometers.  
Width: 112mm. Inner diameter: 12mm  
Outer diameter: 45mm

Part no.	Unit	Price
8500210	1 roll/25m	£3.72

ABR

IMP



### Thermal paper LPQ80, 80m

For **GSI** TympStar™ Pro.  
Width: 80mm  
Inner diameter: 13mm  
Outer diameter: 80mm

Part no.	Unit	Price
8509893	1 roll/80m	£4.23

IMP



### Adhesive thermal paper LPQ80

For **GSI** TympStar™ Pro.  
Width: 82mm. Inner diameter: 25mm  
Outer diameter: 88mm

Part no.	Unit	Price
8509894	1 roll/195 labels	£8.29

IMP



### For Sanibel PRT MPT-II printer, 10m

Part no.	Unit	Price
8013009	1 roll/10m	£3.72



### Label paper roll HM-E200

Part no.	Unit	Price
8519836	1 roll	£3.72







### Otoscope specula

Disposable ear speculums for Heine/KaWe and Welch Allyn otoscopes size 2.5 and 4mm.

Part no.	Description	Unit	Price
8514367	4mm Speculum Economy for Heine/KaWe	100 pcs	£3.72
8514368	2.5mm Speculum Economy for Heine/KaWe	100 pcs	£3.72
8514369	4mm Speculum for Heine/KaWe	100 pcs	£4.65
8514370	2.5mm Speculum for Heine/KaWe	100 pcs	£4.65
8507618	4mm Speculum Genuine for Heine/KaWe	50 pcs	£2.54
8507594	2.5mm Speculum Genuine for Heine/KaWe	50 pcs	£2.54
8108252	4mm Speculum for Welch Allyn	250 pcs	£14.46
8107750	2.5mm Speculum for Welch Allyn	250 pcs	£14.46

AUD BAL IMP LATEX



### Luxamed Otoscope Auris LED 2.5V

Otoscope from German Luxamed. LED light to ensure longevity and great light source, dimable. Fits Heine/KaWe compatible speculums.

Part no.	Unit	Price
8519433	1 pc	£67.64

AUD IMP

### Luxamed Otoscope Auris LED 3.7V USB/PC Charged

Otoscope from German Luxamed. LED light to ensure longevity and great light source. USB charging. Fits Heine/KaWe compatible speculums.

Part no.	Unit	Price
8519434	1 pc	£158.11




### Disposable speculas for video otoscope

Fits the AURICAL™ OTOCam 300, **Interacoustics**® old VIOT and XION™ video otoscopes.

Part no.	Unit	Price
8500470	Baq/50 pcs	£11.67

AUD BAL IMP LATEX



### Video Oscopes specula

Fits the new **Interacoustics**® Viot and OriVision otoscope.

Part no.	Description	Unit	Price
8515083	2.5mm speculum	Baq/100 pcs	£7.44
8511728	4mm speculum	Baq/100 pcs	£7.44

AUD BAL IMP LATEX



### HEINE® mini-c ear light

Part no.	Unit	Price
8530890	1 pc	£67.64



### HEINE® Mini3000® Otoscope

Classic pocket otoscope from Heine in Germany.

Part no.	Unit	Price
8108833	1 pc	£120.91

AUD IMP



### Accessories for Otoscopy

Replacement bulb for Heine mini3000 Otoscope.

Part no.	Description	Price
8515448	Bulb for otoscop Heine mini3000	£24.18
8530892	Bulb for Heine mini-c earlight	£30.44

AUD IMP



### Accessories for Otoscopy

Stylish wall dispenser for ear funnels. Can contain 6 different sizes.

Part no.	Description	Price
8515678	Specula wall dispenser	£39.23

AUD IMP





### EARNimals®

Cartoon faces that are placed on the otoscope. Make ear examinations less scary for kids.

Part no.	Unit	Price
8011501	Baq/5 pcs (assorted)	£8.37

AUD IMP





Battery adaptor

For **Interacoustics** Affinity™, Callisto™, Equinox™, **MedRx** Avant Hit+ and **Siemens** UNITY 2.

Part no.	Description	Unit	Price
8011198	BAA675	1 pc	£37.20
8011199	BAA13	1 pc	£37.20
8011200	BAA312	1 pc	£37.20
8011201	BAA10	1 pc	£37.20
8011202	BAA5	1 pc	£37.20

HIF



Aidapter Interacoustics

Part no.	Description	Unit	Price
8525568	Clear Silicone acoustic couplers	3 pcs	£27.82



Ear coupler sealing wax

Use to create a seal between couplerand hearing instrument.

Part no.	Unit	Price
8500080	1 pc	£7.44

HIF



Hearing Instrument batteries Rayovac

Fits all types of Hearing Instruments.

Part no.	Description	Unit	Price
8520134	Rayovac Extra Zap 10	6 pcs	£3.64
8520135	Rayovac Extra Zap 13	6 pcs	£3.64
8520136	Rayovac Extra Zap 312	6 pcs	£3.64
8520137	Rayovac Extra Zap 675	6 pcs	£3.64

HIF



Battery accessories

Fits all types of Hearing Instruments.

Part no.	Description	Unit	Price
8515452	Battery tester	1 pc	£5.07
8520138	Battery "carrier"	1 pc	£1.61

HIF



Alkaline battery 1.5V AA

Part no.	Unit	Price
8011226	1 pc	£0.68

ABR

AUD

OAE

[supplieshub@interacoustics.com](mailto:supplieshub@interacoustics.com)



Cadzow Industrial Estate  
Low Waters Road, Hamilton  
Lanarkshire, ML3 7QE

**Supplies Hub**  
01698 208205  
supplieshub@interacoustics.com

**interacoustics.co.uk**



**Disclaimer** For advice on compatibility and instructions of use, please consult your operation manual or local distributor. Interacoustics UK assumes no legal liability or responsibility for the accuracy, typing errors or mistakes, correctness or completeness of any information in this brochure. The information in this flyer was correct to the best of our knowledge at the time of printing. Interacoustics UK reserves the right to make changes at any time in specification of products described in this flyer without notice or obligation. MADSEN AccuScreen® and MADSEN OTOflex® are registered trademarks of Otometrics A/S. Echo-Screen®, AuDx®, Navigator® Pro, and ABaer® are registered trademarks of Natus Medical Incorporated. Sanibel™ is a registered trademark of Sanibel Supply®. All other Instrument brands or model names are copyrighted or trademarked by their respective owners and are used only for reference.